

Karma Wheelchair

Datasheet: 00140



Swing-away
Footplates



Carer handbrakes and
Folding backrest



Push-to-lock
Brakes

Karma Ergo Lite Deluxe - Self Propelled

Basic Facts

- Weighing only 10.2 KG, this wheelchair is one of the lightest available
- Equipped with swing away legrests, carer handbrakes, and a folding backrest making it both user and carer friendly

Specifications

- Sizes: 41 cm or 46 cm Wide
- Armrest height: 21.5 cm High
- Seat height: 48 cm
- Backrest height: 42 cm
- Total height: 91 cm
- Weight: 10.2 KG
- Max User Weight: 100 KG

Applications

This wheelchair is commonly used in:-

- Nursing Homes
- Hospitals
- Aged Care Facilities
- Residential

Features

- Folding Backrest
- Lightweight aluminum folding frame
- Fixed seat height
- Ergonomic padded armrest
- Push to lock brakes
- Swing away/removable legrests



P: 5144 4955 F: 5144 6044

E: sales@accessrehabequip.com.au

W: www.accessrehabequip.com.au

156 York Street, Sale, Vic

'Providing Solutions for Life'