

# Backrests



The Mantaray's unique fully adjustable wing tips can be width and height adjusted to give excellent lateral stability.

Offers outstanding **TRUNK AND SPINE STABILITY** to improve performance

Wrap upholstery section for **CONFORMING COMFORT**



Part of the Spex family so is fully compatible with **SPEX ACCESSORIES**

**LIGHTWEIGHT RIGID** aluminium backrest with quick-release hardware

In-built fixed **LATERAL STABILITY**



Laterals are adjustable in:

- Height
- Width
- Angle



## Clinical

Designed for the individual who requires:

- A solid/rigid backrest.
- Ability to self-propel.
- Minimal to moderate postural support.
- Improved spinal and pelvic positioning whilst maintaining propelling efficiency.
- Has mild torso asymmetries.
- Appropriate for clients who have changing postural needs over time.
- Weight capacity with standard shell: 10-15": 90kg/200lb  
16-20": 112kg/250lb
- Weight capacity with heavy duty shell: 10-15": 160kg/350lb  
16-20": 180kg/400lb

## Functional Benefits

- Offers exceptional trunk and thoracic stability.
- Lightweight rigid aluminium backrest with quick-release hardware.
- Minimal contouring, and includes positioning kit for further postural enhancements.
- Modular componentry allows you to 'mix and match' with different lateral options and positioning supports if more forceful support is required.

## Accessories

- Spare outer cover
- Spare incontinence cover

W/C Width	Height	Backrest Code
12" Kidz™	16"	1206-1216-000
14"	18"	1206-1418-000
16"	20"	1206-1620-000
17"	20"	1206-1720-000
18"	20"	1206-1820-000
20"	20"	1206-2020-000