

Tilt- Recline Wheelchair

Datasheet: 00182

Specifications

Code: M4102

Seat Depth:	Seat Width:	Back height:	MUW:
38-50 mm	37-51 mm	50-65 mm	140 KG

Features

- Recline angle 0° - 40°
- Tilt angle 2° - 20°
- Self propel or transit option available
- Removable headrest with adjustable side supports
- Available with self-propelled 24" wheels or transit style 16"
- Carer operated drum brakes standard on transit model and optional on self-propelled
- Alloy construction with water resistant upholstery
- Elevating legrests
- Adjustable armrests
- Adjustable seat depth
- Swing away legrests
- High density memory foam



Applications

The Relax Tilt & Recline Wheelchair is ideal for hospital, nursing home and residential facility environments. Particularly suited to clients that are sitting for long periods of time and may be susceptible to pressure issues.

Description

The Relax is a multi-function tilt and recline wheelchair (including anterior tilt) with a wide range of features and adjustments enabling it to be tailored to individual client needs. The dartex upholstery combined with high density seat foams, ensures a high level of comfort and pressure relief, and postural support.



P: 5144 4955 F: 5144 6044
E: sales@accessrehabequip.com.au
W: www.accessrehabequip.com.au
156 York Street, Sale, Vic

'Providing Solutions for Life'