

Signature Lift Chair

Datasheet: 0283

Signature Lift Chairs are tailored with premium upholstery and durable foams to provide years of postural support, comfort and safety. The Space Saver feature allows the chair to be positioned closer to walls, thereby utilising the best use of your living space. Available in either single or dual motor, the Signature provides support for standing and reclining in comfort.

Features

- Available in Space Saver (Single motor) or Dual Action (Dual motor)
- Depth adjustable seat
- Multi-adjustable padded backrest
- Easy to use backlit handset with USB charger point
- Accessible storage pocket
- Easy clean fabric with resistance to flame, mould, mildew and pilling
- Mink or Heritage

Applications

- Nursing Homes & Aged Care (Low & High Care)
- Residential Homes

Clinical Considerations

FALLS PREVENTION

- Raised posture maximises entry or exit access with mobility aids
- Handset cable runs through side wall pocket to remain in reach even if dropped

OEDEMA MANAGEMENT

- Raised leg positioning to assist with lower limb oedema

CHANGING NEEDS MANAGEMENT

- Adjustable seat depth and backrest upholstery with easy to clean fabric
- Lift feature facilitates reduced assistance as mobility improves in line with clinical rehab goals

DAY BED POSITIONING

- Lie flat posture for improved patient comfort whilst maintaining mobility and independence



Specifications

SPACE SAVER

	Small	Medium	Large
Seat Height	470mm	480mm	510mm
Seat Depth	450mm	490mm	510mm
Seat Width	450mm	480mm	500mm
Seat to Top of Back	590mm	620mm	670mm
Overall Height	1010mm	1050mm	1110mm
Overall Depth	690mm	715mm	740mm
Overall Width	780mm	780mm	780mm
Distance from Wall	140mm	150mm	160mm
User Height Guide	5'2" and under	5'3" to 5'8"	5'9" to 6'2"
Safe Working Load	145kg	145kg	145kg

DUAL ACTION

	Small	Medium	Large
Seat Height	470mm	480mm	510mm
Seat Depth	450mm	490mm	510mm
Seat Width	450mm	480mm	500mm
Seat to Top of Back	590mm	620mm	670mm
Overall Height	1010mm	1050mm	1110mm
Overall Depth	710mm	735mm	760mm
Overall Width	780mm	780mm	780mm
Distance from Wall	430mm	455mm	480mm
User Height Guide	5'2" and under	5'3" to 5'8"	5'9" to 6'2"
Safe Working Load	155kg	155kg	155kg


- * User Height Guide is a reference tool and may not apply to all users
- * Chair measurements may vary +/- 20mm due to upholstery tolerances
- * Seat heights shown are without castors on chair frames
- * Seat width is taken between the thickest sections of armrest padding



Lagoon Suede



Mink Suede



Heritage Vinyl