



# QUASAR

## WHEELCHAIR

Product Code 1595

**Titanium  
Frame**

## A stunning performer on the ground!

The new ultra-lightweight active wheelchair that challenges the wheelchair world with its titanium construction, dynamics and agility.



- Ultra lightweight wheelchair.
- Titanium rigid frame.
- Frame weight: 4.4 kg without rear wheels (*7kg with wheels*).
- Available seat width: from 34 to 44 cm.
- Titanium open foot-rest.
- Aluminium clothes guards with carbon fibre wings.
- Foldaway brakes.
- New 24-spoke Ultralight wheels.
- Available in fixed or adjustable rear axle



New side guard bracket includes backrest fold down mechanism and seat to back angle adjustment.



Adjustable height front forks with wheel size options.



Folds easily and compact for storage or transportation.

### Fixed Axle



### Adjustable Axle



Adds 400gm



**Accessory:** Max Grip Ergo push rims for ultimate push-power and comfort for the users hands.



Unique single-sided anti-tipper. Easily swings away when not needed.

[www.accessrehabequip.com.au](http://www.accessrehabequip.com.au)

Call (03) 5144 4955

Fax (03) 5144 6044

[sales@accessrehabequip.com.au](mailto:sales@accessrehabequip.com.au)

156 York St, Sale VIC 3850

PO Box 1376, Sale VIC 3850

# QUASAR

## ORDER FORM

(Please copy for multiple use)

1595

DATE: / /

REQUEST FROM:

ORGANISATION:

POSTAL ADDRESS:

PHONE:

FAX:

PURCHASE ORDER NUMBER / FUNDING DETAILS:

DELIVER TO:

ORGANISATION:

PHYSICAL ADDRESS:

ATTENTION:

FOR (USERS NAME):

M / F

☐ **QUASAR FIX**

(Completely fixed)

☐ **QUASAR ADJ**

(Adjustable rear axle)

**Key**

- ☐ standard  
☐ option

**L** Seat width ☐34 ☐36 ☐38 ☐40 ☐42 ☐44  
**P** Seat depth ☐36 ☐38 ☐40 ☐42  
**HS** Back-rest height ☐27 ☐29 ☐31 ☐33 ☐35 ☐37  
**HA** Front seat/floor height ☐46 ☐48 ☐50  
**HP** Rear seat/floor height ☐40 ☐42 ☐44 ☐46 ☐48  
**DP** Foot-rest/seat distance ☐37 ☐39 ☐41 ☐43 ☐45 ☐47 ☐49

cm

**Frame colour**

- ☐ Shiny Titanium ☐ Bronzed ☐ Pearl white  
☐ Glossy black ☐ Anthracite ☐ Metallic grey ☐ Cyclamen Light Purple

☐ ..... (see "Colours and Accessories" leaflet)


①

**Back-rest upholstery**

- ☐ Breathable  
☐ Nylon (black colour)

②

# QUASAR




③  
**Single tubular open foot-rest**



☐ Closed

☐ **Without pushing handles**



☐ **Fixed pushing handles**

④

⑤  
**Clothes guards**


☐ Aluminium straight



☐ Aluminium with carbon fibre wing

⑥  
**Brakes**

☐ Standard



☐ Foldaway



**QUASAR FIX**


**7 kg weight Set-up**

- Wheelchair 38x40x31 cm
- Nylon upholstery
- Shiny titanium frame colour
- Open foot-rest
- Without pushing handles
- Aluminium clothes guards with carbon fibre wing
- Foldaway brakes
- Ultralight wheels
- Pneumatic high pressure smooth tyres
- Ultralight push rim
- 80 mm front castors

Max load: 100 kg

⑦  
**Distance between the hub centre and the back-rest**

- ☐ 7 cm  
☐ 9 cm (standard)  
☐ 11 cm  
☐ 13 cm  
☐ 15 cm



⑧  
**Back-rest inclination**


93° 89°  
91°



- ☐ 91° (standard)  
☐ 93° ☐ 89°

⑨  
**Revolving anti-tip device**

☐ RIGHT or ☐ LEFT



**Rear wheel camber**

- ☐ 1° ☐ 4°

**Rear wheels**

- ☐ 24" Ultralight (24 spokes) ☐ 24" (36 spokes)  
☐ 24" Light (24 spokes)  
☐ 24" Star 6x3 (18 spokes) ☐ 24" off-road  
☐ 24" Spinerger (18 raggi)  
☐ 24" black Racer  
☐ 24" grey Racer  
☐ 24" red Racer  
☐ 24" blue Racer



**Tyres**

- ☐ Pneum. high press. smooth ☐ Pneum. Marathon Plus  
☐ Pneum. high press. plugged ☐ Solid

**Push rims**

- ☐ Aluminium Ultralight (\*\*)  
☐ Aluminium  
☐ Chromium-plated steel  
☐ Rubber coated  
☐ Steering (\*)  
☐ Tween Grepp  
☐ Max Grepp  
☐ Max Grepp Ergo(\*)  
☐ Silicon covers

(\*) : not available with Racer wheels  
(\*\*) : available with Ultralight wheels only

**Spoke guards**



- ☐ Transparent  
☐ White  
☐ Black  
☐ Multicolour  
☐ Black plumed red  
☐ Black plumed blue  
☐ black plumed green

not available with Racer wheels

**Front castors**

- ☐ 80 (3") solid with lights  
☐ 80 (3") solid  
☐ 100 (4") solid with lights  
☐ 100 (4") solid  
☐ 100 (4") solid with alum. rim  
☐ 125 (5") solid  
☐ 125 (5") solid with alum. rim  
☐ 150 (6") pneumatic  
☐ 150 (6") solid

**Pockets**



- ☐ Small under the seat ☐ Big under the seat ☐ on the back-rest

⑭

**Coloured set**

(Rim, hub, forks, clothes guards support)

- ☐ Shiny red  
☐ Shiny blue  
☐ Bronzed

⑮  
**Lined cushion**  
(foam polyurethane)

- ☐ 3 cm high  
☐ 5 cm high

**Other options**

- 17 ☐ Male velcro strip on seat  
18 ☐ Female velcro strip on seat  
19 ☐ Holdall