


Quadform Mattress

Datasheet: 00138



Fully Supportive Immersion



Four layers of Pressure-relieving foams



Vapour Permeable, 4-way stretch cover

Basic Facts

- 4 Layers of quality foam specifically designed to provide unparalleled levels of comfort and pressure care
- High Profile (18 cm) thickness
- Firm side bolsters to assist with transfers
- Flexible sections for use with profiling beds

Commonly Available Sizes

- Single 90 cm x 200 cm
- King Single 107 cm x 200 cm
- Double* 137 cm x 200 cm
- Queen* 150 cm x 200 cm
- *Custom sizes available on request

Description

- Constructed from 4 layers of Dunlop foams in graduated density for optimum support
- Covered in 4-way stretch Carrflex Vapour permeable pressure care material
- The top layer of foam is plush memory foam for ultimate comfort and pain relief

Applications

This mattress is commonly used in:-

- Nursing Homes
- Hospitals
- Aged Care (Low Care and High Care)
- Residential
- Palliative Care

Features

- Polyurethane cover
- Cover has a zipper for ease removal
- Cover is 100% waterproof, aiding hygiene and cleaning processes



P: 5144 4955 **F:** 5144 6044

E: sales@accessrehabequip.com.au

W: www.accessrehabequip.com.au

156 York Street, Sale, Vic

'Providing Solutions for Life'