

Regency Care-Chair

Datasheet: 00137



'The Chair that Cares'

Fully moveable sides and tilt-recline features make transferring and moving as easy as 1,2,3!



Description

- The Regency Care Chair is the market leader in providing ultimate comfort and pressure relief for the user, while helping to make the caring role easier. Cleverly designed features such as the patient lifter friendly castor base provide endless benefits to all involved.

Features

- Forward Transfer friendly with fold-away legrest and footplate
- 175mm Manoeuvrable castors for ease of movement
- Profiled head cushion for increased comfort
- Side-transfer friendly with patient lifter friendly castor base
- Superior range of tilt and 'Virtual Flat' recline features
- 2 position armrests
- Self-adjusting backrest
- Multi-positional Seating configuration
- Armrests can be fitted with a Meal tray
- Large range of options
- Larger models available

Model: R200D pictured

Applications

This Chair is commonly used in:-

- Nursing Homes & Hospitals
- Aged Care (Low Care and High Care)
- Residential
- Palliative Care

Specifications

- | | |
|------------------------------------|---------|
| • Armrest height (from seat) | 175 mm |
| • Backrest height (from seat) | 680 mm |
| • Footrest - Extendable | + 85 |
| • Overall height (backrest up) | 1260 mm |
| • Overall length (supine position) | 1800 mm |
| • Overall width | 660 mm |
| • Seat depth | 520 mm |
| • Seat height (normal position) | 630 mm |
| • Seat width between arms | 460 mm |
| • Max User Weight | 180 KG |



P: 5144 4955 F: 5144 6044
E: sales@accessrehabequip.com.au
W: www.accessrehabequip.com.au
156 York Street, Sale, Vic

'Providing Solutions for Life'