

Posture Plus Chair

Datasheet: 00143



Matching upholstered arms



Fully ergonomic back support



Rear wheels make repositioning easy

The Posture Plus Chair

Description

The Posture Plus is the final word in patient and day room seating. Featuring an ergonomically sculptured seat shell and seat foams, as well as integrated head support, the Posture Plus offers outstanding user comfort, even when used for extended periods.

Features

- Ergonomically sculpted seat shell and seat foams
- Height adjustable legs
- Available in a wide range of vinyls
- Arms can also be upholstered to match
- Comes standard with transit wheels on the rear legs for ease of repositioning
- Comfortable and supportive high back
- Black, silver and beige frame powdercoat available as standard

Applications

This Chair is commonly used in:-

- Nursing Homes & Hospitals
- Aged Care (Low Care and High Care)
- Residential Homes

Specifications

- | | |
|-------------------------------|--------------|
| • Armrest height (from seat) | 175 mm |
| • Backrest height (from seat) | 680 mm |
| • Overall width | 640 mm |
| • Overall height | 1020-1135 mm |
| • Overall depth | 740 mm |
| • Seat height | 420-535 mm |
| • Seat width | 540 mm |
| • Seat depth | 470 mm |
| • Max User Weight | 180 KG |



P: 5144 4955 F: 5144 6044
E: sales@accessrehabequip.com.au
W: www.accessrehabequip.com.au
156 York Street, Sale, Vic

'Providing Solutions for Life'