

# Thompson Chair

Datasheet: 00144



Padded arms can be upholstered to match



Fully Ergonomic back & seat support



Height adjustable legs

## *The Thompson Utility Chair*

### Description

The Posture Plus is the final word in patient and day room seating. Featuring an ergonomically sculptured seat shell and seat foams, as well as integrated head support, the Posture Plus offers outstanding user comfort, even when used for extended periods.

### Features

- Ergonomically sculpted seat shell and seat foams
- Height adjustable legs
- Available in a wide range of vinyls
- Arms can also be upholstered to match
- 25 mm Steel tubular frame
- Contoured and supportive back and waterfall seat front to reduce pressure behind the knees
- Black, silver and beige frame powdercoat available as standard

### Applications

This Chair is commonly used in:-

- Nursing Homes & Hospitals
- Aged Care (Low Care and High Care)
- Residential Homes
- Waiting rooms and Showrooms
- Clinic Rooms

### Specifications

- |                               |              |
|-------------------------------|--------------|
| • Armrest height (from seat)  | 175 mm       |
| • Backrest height (from seat) | 680 mm       |
| • Overall width               | 565 mm       |
| • Overall height              | 825 - 950 mm |
| • Overall depth               | 600 mm       |
| • Seat height                 | 460 - 585 mm |
| • Seat width                  | 540 mm       |
| • Seat depth                  | 480 mm       |
| • Max User Weight             | 120 KG       |



**P: 5144 4955 F: 5144 6044**  
**E: [sales@accessrehabequip.com.au](mailto:sales@accessrehabequip.com.au)**  
**W: [www.accessrehabequip.com.au](http://www.accessrehabequip.com.au)**  
**156 York Street, Sale, Vic**

*'Providing Solutions for Life'*