

Evolution Chair

Datasheet: 00145



Elevating
legrest



Drop-sides for
easy transfers

Low-friction Castors

'The Patient and Carer-friendly Care-chair'

Description

- The Evolution Supreme Chair has been designed as a low-maintenance, user friendly Care-chair with a large range of seating adjustments, allowing it to suit almost any patient. Soft touch levers for ease of operation, a lower seat height, and low friction castors help to care for not only the patient, but the care-giver as well.

Features

- Covered in Carflex multi stretch cover to the contact areas
- Low seat height for easier stand assist transfers
- Soft-touch levers for easier operation of tilt, recline and the elevating legrest
- Drop-sides standard for ease of side transfers
- Low friction 5" castors with directional lock
- Curved handles at rear to assist with ease of moving and repositioning
- Multi-positional Seating configuration
- Visco-elastic foams distribute weight over a larger area without compromising postural support or the natural spine alignment
- Accessories such as Meal Tray, Foot rests and Legrest extensions available
- Forrest or Coffee (pictured) colours available

Applications

This Chair is commonly used in:-

- Nursing Homes & Hospitals
- Aged Care (Low Care and High Care)
- Residential
- Palliative Care

Specifications

	Small	Medium	Large
Seat Depth	485 mm	585 mm	585 mm
Back Height	800 mm	800 mm	800 mm
Seat Width	455 mm	485 mm	585 mm
Seat to Floor	620 mm	620 mm	620 mm
Overall Width	700 mm	730 mm	830 mm
MUW	145 KG	145 KG	190 KG



P: 5144 4955 **F:** 5144 6044
E: sales@accessrehabequip.com.au
W: www.accessrehabequip.com.au
 156 York Street, Sale, Vic

'Providing Solutions for Life'