


Shoprider

Rocky 4

ASK FOR THE BRAND **SHOPRIDER**



ROCKY 4

SHOPRIDER®
mobility
www.shoprider.com.au

RELIABLE MOBILITY GUARANTEED

Shoprider

Rocky 4

SHOPRIDER®
www.shoprider.com.au mobility

The SHOPRIDER™ **Rocky 4** is an outstanding model from the extensive range of mobility scooters available from SHOPRIDER. The **Rocky 4** must be seen to be believed, it redefines the concept of reliable mobility. The **Rocky 4** comes with 330mm wheels, Airflow Comfort Seat and all the features you expect from the SHOPRIDER range, fully rotating lockable seat, suspension, adjustable steering tiller and multi position lift up armrests.

The reliability of this impressive scooter by SHOPRIDER™ is further enhanced by the rear wheel drive and electromagnetic brake.

SPECIFICATIONS

ROCKY 4

Length	1430mm
Width	680mm
Height	1300mm
*Safe Climbing Angle	9 Degrees
*Max Climbing Angle	12 Degrees
*Maximum Load <u>Single user</u>	225Kgs
Charger	8Amp
Motor Size	2.0Hp
Maximum Speed	15Kph
Battery capacity min	2 x 60 Amp
Maximum range#	up to 50km
Wheel Size	330mm
Turning Radius	1800 mm
Weight (with battery)	146kg

FEATURES:


Airflow Comfort Seat
 Loop Steering Tiller
 Adjustable Steering Tiller
 Adjustable width arm rests
 Sliding Rotating Seat
 Variable Seat Height
 Battery Power indicator
 Double Head Lights
 Rear Lights & Indicators
 Hazard Lights
 Front & Rear Bump Pads
 High & Low Speed
 Shopping Basket
Comprehensive Warranty

WARNING: Do not exceed 10kph on public roads or paths

IMPORTANT: Always be sure to read your owner's manual thoroughly before operating your machine.

Driving range varies with the factors shown below

- (1) The weight of the occupant
- (2) Level ground & Ambient temperature
- (3) Battery condition / Type of charger
- (4) Driving Style



*(Single user)

Shoprider is listed and complies with the regulations of the Therapeutic Goods Association of Australia as well as the regulatory bodies and standards in America (FDA) and Europe (CE)

SHOPRIDER®
www.shoprider.com.au mobility

Your Authorised dealer

For your nearest Distributor call 1800 672 633
 Or visit our Web Site at www.shoprider.com.au