

Meuris Recliner Club Chair



Features & Benefits

- ◆ Designed with therapists and care providers for ultimate patient comfort.
- ◆ Vehicle impact tested to Australian Standards and Nationally approved for vehicle transport.
- ◆ Vehicle restraints which conform to Australian Standards.
- ◆ Gas assisted seat tilt, leg raise and back recline.
- ◆ Removable arm rest for side transfers and handling.
- ◆ 4 way stretch fabric for minimal shear.
- ◆ Front wheel steering.
- ◆ Superior ergonomics and comfort through the use of Pressure Mapping of the Meuris adjustable air cell system.
- ◆ Removable tray table and dual footplates included.
- ◆ 6" lockable castor wheels and ergonomic handles for easier portability and manoeuvring.
- ◆ Handles also provide for ease of recline and increased user comfort.



Specifications

Back Width:	48cm (19")
Seat-to-Floor:	62cm (24.5")
Seat Depth:	50cm (20")
Seat Width:	48cm (19")
Top of Back to Floor:	120cm (47.25")
Weight Capacity:	125kg

