


6001 Series Bed

Datasheet: 00147



Description

The 6001 Series Bed is the standard base model electric bed, designed for low care patients. Standard care features have been included, such as electric backrest, electric knee-break, and a 200KG Safe working load.

By using 125 mm braking and directional castors, confusing directional lock systems have been removed, and hence the bed is simple and easy to manoeuvre.

Features

- Electric Hi/Low
- Electric Backrest
- Electric Kneebend
- SWL 200KG
- 125 mm Braking & Directional Castors
- Minimum Height: 270 mm
- Maximum Height: 700 mm (add 300 mm to above heights for Central locking)
- Overall width: 995 mm
- Overall Length: 2060 mm
- Sleeping surface: 1990 x 870 mm

Applications

The 6001 Series Bed is commonly used in:-

- Nursing Homes
- Hospitals
- Aged Care (Low Care and High Care)
- Residential
- Palliative Care

Accessories and Options

- Side rails (vertical or horizontal, or clamp on)
- Padded side rail protectors
- Self-help or I.V Poles
- Battery back-up
- Attendant control
- 125 mm Central locking castors
- CPR Emergency Release
- Head and Footboard style and finish options (Custom options available on request)



P: 5144 4955 F: 5144 6044

E: sales@accessrehabequip.com.au

W: www.accessrehabequip.com.au

156 York Street, Sale, Vic

'Providing Solutions for Life'