

Trust Outdoor Walker

Datasheet: 00151


Height adjustable handles and dual action brakes

Large wheels, perfect for outdoors



Faux leather seat

Sprung wheels for increased comfort.



Folds, and stands alone when folded.

Description

This sleek and lightweight walker, with its large wheels is the practical mobility solution for outdoors. The handles are height adjustable, and dual locking for versatility. The spring system on the wheels provides a smooth and comfortable ride. It can be folded, and is able to stand alone once folded, making it easy to move in and out of the car.

Features

- 360° Front wheel provides maximum manoeuvrability.
- Dual action brakes
- 'X' Style frame folds, and stands alone when folded.
- Sprung wheels for comfortable pushing - Reduces jarring.
- Height adjustable, dual action handles.
- Faux leather seat.
- Available in beige or black frame.

Applications

This walker is most suitable for outdoor use, but can also be used indoors.

Specifications

Handle Height:	78 - 96 cm
Width between handles:	41 cm
Seat Height:	59 cm
Length:	64 cm
Front wheel Dia:	24.5 cm
Rear wheel Dia:	20 cm
Overall Width:	57 cm
Weight:	6.5 KG
Maximum User Weight:	130 KG



P: 5144 4955 **F:** 5144 6044
E: sales@accessrehaequip.com.au
W: www.accessrehaequip.com.au
156 York Street, Sale, Vic

'Providing Solutions for Life'