

E-VIP Commode

Datasheet: 00152



Description

- The Ocean E-VIP Commode is a revolution for patient showering requirements. It is an electrically height adjustable, multi-purpose shower commode designed specifically for reducing the strain carers often experience when assisting patients while showering. The E-VIP has a unique tilt function to assist carers with showering even the most passive of patients whilst also working towards the prevention of back injury.

Features & Benefits

- Exceptional range of tilt and height adjustments achieved by simply pressing a button.
- Seat can either be tilted anteriorly or posteriorly (-5° to 35°) to enable greater access to patients with decreased upper body control.
- All electric components are sealed for safety and durability.
- Chest and pelvic belts available.
- Optional wider XL armrests to increase the seat width.
- One rear castor is equipped with directional lock for ease of transporting.
- Adjustable headrest can be set horizontally or vertically to users requirements.

Applications

- The Ocean E-VIP Commode is designed specifically to:-
- Reduce carer strain.
 - Provide a safer and more comfortable showering experience for the carer and the patient.

Specifications

Height Range:	500-1050mm
Seat Width:	480mm
Unit Weight	39 KG
MUW	150 KG



P: 5144 4955 F: 5144 6044
E: sales@accessrehabequip.com.au
W: www.accessrehabequip.com.au
156 York Street, Sale, Vic

'Providing Solutions for Life'