

2001MKII Series Bed

Datasheet: 00148




Wall Protector Buffer



Sheet Metal Mattress Base



Braking & Directional Castor



Description

The 2000 Wide Series Bed is designed with the larger patient in mind. With a 220KG SWL and 1070 mm width, it provides for a wide variety of patient needs.

Features

- Electric Hi/Low
- Electric Backrest
- Electric Kneebend
- SWL 220KG
- 125 mm Braking & Directional Castors
- Minimum Height: 370 mm
- Maximum Height: 700 mm (add 300 mm to above heights for Central locking)
- Sleeping surface length: 2010 mm
- Manufactured in Australia.

Applications

The 2001 MKII Wide Series Bed is commonly used in:-

- Nursing Homes
- Hospitals
- Aged Care (Low Care and High Care)
- Palliative Care
- Residential

Accessories and Options

- Side rails (vertical or horizontal, or clamp on)
- Padded side rail protectors
- Self-help or I.V Poles
- Wall buffer-bar
- Battery back-up
- Attendant control
- 125 mm Central locking castors
- CPR Emergency Release
- Head and Footboard style and finish options (Custom options available on request)



P: 5144 4955 F: 5144 6044

E: sales@accessrehabequip.com.au

W: www.accessrehabequip.com.au

156 York Street, Sale, Vic

'Providing Solutions for Life'