

BREEZY[®]

BasiX²



*Improved
adjustability
&
rigidity combined
with a modern
look!*



**Provides an even broader
range of clinical solutions**

 Breezy



BasiX²

Key Product Features:

Minimize Inventory:

- **1 Chair - Multiple Depths**
Adjust Seat Depth by sliding cross brace 16" to 18" deep
- **1 Chair - Enhanced wheel adjustability to suit more clients**
4 vertical positions, 2 active and 2 passive horizontal wheel positions
- **Self Propel or Transit wheels available**
- **Tension adjustable backrest now extendable from 41 to 46 cm**

Upgraded Features:



Axle Plate



Caster Housing



Wheel Lock



Upholstery

New Options:



Height Adjustable arm rests



Adjustable Handles

Technical Specifications:

Max. User Weight:	125 kg	Back:	7° bend (std), ½ folding
Chair Weight:	14.9 kg	COG:	+3.5; - 5 cm
Transport Weight:	8.4 kg	Armrest: (std)	20 - 30 cm (covers short & full length pad in one)
Overall Width:	SW + 19 cm / SW + 20 cm (drumbrake version)	Armrest height adj.:	20.5 - 23 cm
Folding Width:	Max. 30 cm	Height Adj. Push Handles	+12 cm
Seat Width:	14, 16, 18 & 20"	Hangers:	70° (lower leg length adjustable from 35 - 50 cm) Angle adjustable footplate (+/- 15°)
Seat Depth:	16 - 18" - adjustable	Colour:	- Robust Hammerite - makes scratches less visible - Red (on selected SKU's)
Seat Height:	48 cm (adj. 42, 5 - 51 cm)	Crash Test Approved	Yes
Back Height:	41-46 cm		
TGA#:	145834		
Warranty:	Limited Lifetime		

Sunrise Medical Pty Ltd
6 Healey Circuit
Huntingwood NSW 2148
Ph: 02 9678 6600
Fax: 02 9831 2244
orders@sunrisemedical.com.au
www.sunrisemedical.com.au